3. FIT

- a) Trim, if required, at the trim sections **C** with scissors while the orthosis is hot.
- b) Fold and stack the last level onto the next level and compress to form a comfortable cuff.

4. APPLY

- a) Rest the patient's elbow on a flat surface with the thumb abducted. Fit the spica over the thumb and align the radial hinge **D** along the radius.
- b) Slide the base of the spica well down around thumb and ensure the thumbhole bridge **I** is snugly positioned.
- c) Stretch & wrap the palmar cuff **E** around the hand and align the closure tab **G** onto the closure fixing zone **F** & press into place.
- d) Stretch & wrap wing **B** around the arm to overlap wing **A** over its full length. Align & secure yellow dot **N** with yellow dot **M**. Ensure semi-circular tab **H** overlaps palmar area **J**.
- e) Easy-fix tabs **L** will temporarily secure wing in place. Confirm the radial hinge **D** is still aligned along the radius.

5. MOLD

- a) Roll the fabric liner **L** protruding from the Spica, back over itself to form a padded edge, typically ensuring patient can flex IP joint.
- b) Wrap an elastic bandage (not included) firmly but not tightly to conform the spica and entire orthosis. A bandage wetted with iced water accelerates cooling.
- c) Manipulate the device into the most appropriate position for healing, typically using Arm-wrestler grip and cupping for dorsal support up the forearm proximally.

6. SET

- a) Rest patient's arm on a table until the orthosis is rigid. Remove the elastic bandage and easy-fix tabs. Verify that the wing **B** and palmar cuff **E** are free to allow for subsequent adjustment.

7. SECURE

- a) Apply the Compliance-Lock strap to the wrist area of the orthosis. Thread the strap through the slot in the liner. Thread it through the D-ring & tighten for a snug fit.
- b) To lock this strap, push the white fastener peg through the nearest matching hole on the strap and snap to lock in place. The strap can be loosened a pre-determined distance e.g. to allow for swelling. To tighten it further break link between fastener and D-ring. To remove the compliance strap, cut the white peg with scissors.
- c) Fit forearm tension strap and tension for a snug fit.

Note: the physician should advise the patient on how to adjust the straps for a snug fit.

Note: After removal from the patient, the orthosis can be reformed with a heat gun in a localised area or fully re-activated in an oven. See www.FastFormMedical.com for details..

8. REHAB

- a) To prepare orthosis for rehab use, remove it & flex the radial hinge **D** fully for easy opening.
- b) Likewise flex the Palmar Cuff Hinge

GETTING IT WET

If allowed by the physician, it is acceptable for the orthosis to become wet during use. It may take a number of hours to become completely dry. If allowed by the physician to loosen or remove the orthosis, the patient may dry it using a hair dryer on a low-heat setting.

CLEANING INSTRUCTIONS: To clean the orthosis, rinse it in warm water. Note: it is not recommended to use household cleaners or personal hygiene products directly on the orthosis.

STORAGE: This device is not sterile and should be stored between 4 - 35°C (39-95°F). It must also be protected from strong direct sunlight.

WARRANTY: FastForm Research Ltd covers defects in workmanship and materials for up to 12 weeks after this device is dispensed. Damage caused by normal wear & tear, abuse or alterations is not covered.

For further details on heating & application method, see:

www.breg.com/Ffastform/HeatingTimes or www.FastFormmedical.com

-Custom mouldable for single patient use. -Not made with natural rubber latex -Rx only

Notice: while every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this device. If a serious incident occurs with this device, please report it to the Manufacturer and your Competent Authority.

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Instructions for Use Thumb Spica Multifunction Orthosis TSMO

Before using the device, please read the following instructions completely and carefully. Correct application and care is vital to proper functioning of the device.

INTENDED USE / INDICATIONS:

- Treatment of non-displaced fractures of scaphoid and bones of the MCP joints which may co-occur with fractures of the wrist.
- Ulnar ligament injury, inflammatory conditions such as de Quervain syndrome.
- Provide support in the post-healing rehab phase.

CONTRAINDICATIONS:

Complex, compound and/or comminuted fractures.

WARNINGS AND PRECAUTIONS:

Device to be fitted by a qualified healthcare professional. Patients should consult their medical professional immediately if they experience any pain, swelling, sensation changes, or any unusual reactions while using this device. This Device should not be used on a patient with polyester or elastane allergy. Ensure the orthosis is not too hot prior to application. Do not apply over open wound.

Instruct the patient on correct adjustment of straps to maintain snug fit.

Remove the device before any re-molding with a heat-gun, to avoid burning the patient.

The orthosis will lose its shape above 57°C/135°F or under extreme sunlight. Cover with stockinet or clothing.

Avoid water or vapor above 40°C/104°F for prolonged periods. Avoid heat sources (e.g. ovens, flames, sunny car windows, very humid conditions) for long periods.

This device reduces mobility & dexterity. Do not operate heavy equipment.

Patients should consult their medical professional before traveling by airplane.

1. SELECT & GAUGE (See matching photos on packaging insert)

- a) Select Thumb Spica size based on sizing chart. Ensure thumb liner **L** extends out of spica.
- b) To pre-size the orthosis, place it on the patient's arm, aligning the distal dorsal end parallel to the knuckles and measure from the end of the orthosis to the crease on the inside of the elbow. It should be approximately 3 finger-widths away from the crease. Note how many strips (if any) need to be cut off after heating.

2. HEAT

Lay the orthosis fabric side down with the thumb section flat on the FastForm rack, ensuring the plastic does not touch the sides. Heat the orthosis until it becomes soft / malleable. If not fully malleable, re-heat in 1-minute increments. DO NOT heat for a total time greater than one hour

Convection ovens (preheated) = 5 -10 minutes @ to 75-85°C (167-185°F)

Nuwave ProPlus/Elite Ovens = 6 mins @ 220°F

Nuwave Primo = 6 mins @ 260°F

Waterbath Not suitable

For more detail see instructions on heating device/or QR Code.

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FastForm Research Ltd., 4 Lear Ave., Airport Business Park, Killowen, Waterford, IRELAND.
Breg, Inc., 2885 Loker Ave. East, Carlsbad, CA 92010, USA